

supplementary File 1

Table S1. List of primers used for qRT-PCR

Gene name	Primer sequence (5'-3')	
miR-302d	F: GCGTAAGTGCTTCCATGTTTGA	RT:CCAGTGAGCAGAGTGACGAGGACTC GAGCTCAAGCTTTTTTTTTTTTTTTTTTTTAC
miR-371-5p	F: GCGACTCAAAGTGTGGGGGC	RT:CCAGTGAGCAGAGTGACGAGGACTC GAGCTCAAGCTTTTTTTTTTTTTTTTTTTAG
miR-372	F: GGCGAAAGTGCTGCGACATTTGA	RT:CCAGTGAGCAGAGTGACGAGGACTC GAGCTCAAGCTTTTTTTTTTTTTTTTTTACA
U48	F: TGATGACCCCAGGTAAGTCTG	RT:CCAGTGAGCAGAGTGACGAGGACTC GAGCTCAAGCTTTTTTTTTTTTTTTTTTGG
Universal R	CCAGTGAGCAGAGTGAC	